InBodyBAND2

Quick Start Guide

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

※ Check our website inbody.com to download the detailed guide.

1. Powering Up and Charging

- 1) Press and hold(5sec) the electrodes to turn on the InBodyBAND2. (You will see InBody in the center of the screen)
- 2) If the product does not turn on, connect it to the charger. If the charging terminal fits well, a dotted line will appear on the screen.


- ※ A full charge requires approximately 1 hour.
- ※ The App will not work while charging.

2. Turning on Screen and Changing Screen Mode

- 1) Turn the wrist inward like looking at a wristwatch, or press and hold the electrodes instantly with two fingers to turn the screen on.
- 2) While the screen is on, tap the screen or press and hold the electrodes instantly with two fingers to change the screen.


- ※ If your fingers are too dry, the screen may not change.
- ※ Press "Setup > Device Settings" at the upper left of the main screen of the InBody App to change the screen mode.

3. Installing and Syncing with the App

- 1) Search "InBody" in the App Store or the Google Play Store and click download.
- 2) Open the app and follow the guidance to sign in.
- 3) Select "Setup > Device Settings > Device > InBodyBAND2" at the upper left of the main screen of the app while Bluetooth is on.
- 4) Step and sleep data can be synchronized through scrolling down the main screen of the app.

4. Performing InBody Test

- 1) Press the electrodes when a finger icon appears on the measuring screen.
- 2) Press and hold the electrodes for a long time on any screen to start performing InBody Test.

- ※ Keep your fingers away from the back of your hand.
 - * You may want to insert a picture to indicate fingers not touching the back of the hand.
- ※ If your fingers and wrists are dry, it may affect the measurement.
- ※ If the firmware version ends with an "N" (ie, V038N), the band can be connected to InBody 970.

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